

## **Cura Alarms - Patient Activated Monitors**



- Patient activated alert monitor assists with falls prevention
- Attached to patients clothing, patient activates by pulling string which detaches magnet and sets off alarm
- Stand alone system is simple to use
- Durable and reliable
- Portable and lightweight
- Supplied with 9 volt battery

## Warranty Information

• 1 Year



Code Display Name

BEA007460 Pull String Monitor



## ndis

**Item Code:** 

05\_222704219\_0123\_1\_2

Personal Alarm - Composite:

05\_427\_0103\_1\_2

Click **here** for more information

## **Compliance**

No applicable compliance standard