

## Pride Power Assist



**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>

The **Pride Power Assist** is designed to aid caregivers in pushing a manual wheelchair in otherwise strenuous situations such as up steep hills, inclines and loose gravel.

- The Power Assist attaches directly to manual wheelchairs easily.
- It features an easy-to-use squeeze trigger mechanism, along with a speed dial, enabling users to customise speed settings.
- The Power Assist is designed to work in conjunction with most wheelchairs with seat widths of 16" – 20".
- Weight capacity of 120 kg / 264 lbs. (wheelchair and user combined).
- The Power Assist has a maximum speed of up to 5 km/h (3.1 mph).

✔ Contracted DVA product



Give your manual wheelchair a heavy-duty boost with the **Pride Mobility Power Assist HD**. It comes with an easy-to-use controller with speed dial, battery indicator and hand-controlled lever.

- The Pride Power Assist is designed to aid caregivers in pushing a manual wheelchair.
- The PowerAssist attaches directly to manual wheelchairs easily.
- Designed to work with manual wheelchairs with widths 406 mm - 610 mm (16" – 24").
- Range of up to 16 km (10 miles) between charges.
- Weight capacity of 182 kg / 400 lbs. (user and wheelchair combined).
- Maximum speed of up to 6.4 km/h (4mph).
- Dual-wheel design for added traction and stability.

\*Pictured mounted to wheelchair (not included)

Code	Display Name	Safe Working Load	Maximum Speed	Range
WCA858501	Power Drive Unit - PowerAssist - Power Drive Unit - Pride	120 kg	5 km/h	16 km
WCA858511	Power Drive Unit - Power Assist/Stroll HD - (As per Script)	182 kg	6.4 km/h	16 km