

magic Mobility

PACK UP AND HIT THE ROAD

The Fold Forward Backrest is a great feature when travelling with your powerchair. A simple quick release mechanism allows you to fold your chair to nearly half the height.

TAILORED FOR YOU

Select the right headrest, backrest, armrests, legrests, seating and cushioning to suit your lifestyle.



TAKE CONTROL

An extensive range of advanced control systems and environmental controls are available. Easy to program and use.

FULL LIGHTING PACKAGES

Been out having too much fun and it's getting late? Put your lights on.











FRONTIER V6 COMPACT 73

Compact 73 provides users with the ultimate spaces. The Frontier family has proven itself for nearly two decades. The C73 handles kerbs and takes potholes in its stride, giving you improved freedom in and around the streets.

Basic specifications:*

- · Length: 955mm
- Width: 655mm
- Seat to floor height: 435mm
- Base turning radius: 530mm
- Weight limit: 182kg (400lbs), 155kg (340lbs) with lift and/or tilt
- Top speed: 10kph (6.2mph)
- Motors: 800W Gear in Line
- Wheels: Hybrid 14x3" black or grey (solid or pneumatic)
- Castors: 200x50mm black or grey (solid or pneumatic)
- Gel Batteries
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Rehab seat pan (from 12-24" wide)
- Crash test compliant tie down points

Rehab positioning:

- Power options: 300mm (12") seat lift, tilt (45 degrees), lift and tilt, transfer lift and tilt only, transfer tilt with lift and tilt, power recline and more
- Seating and backrests: MPS, rehab backrest, adjustable seat pans and seating
- Headrests and mounts: MPS, rehab
- Legrests: center post, swingaway, power or gas and manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any specific requirement
- A full range of armpads, straps, pads, harness, restraints, supports are available

The right wheelchair for:**

- Every day off-road recreation
- Urban outdoor living and working
- Roomy indoor spaces
- Tight indoor spaces

Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Retractable docking pin
- 24V LED lighting kit
- 12V accessory socket and 2 x USB charger
- Wheelchair jack
- Tray tables
- Luggage rack, bag hooks, bags, pockets, keep safe secret boxes
- Special controls: attendant. head, chin, foot, sip and puff. environmental controls
- RAM® mounting system (tablet and swing away tablet holder, mobile phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder)

Base options:




SUNRISE YELLOW BASE GREY HYBRID TYRES




FOREST GREEN BASE GREY HYBRID TYRES



SPACE BLUE BASE BLACK HYBRID TYRES



not be driven through water, into rivers, creeks or the sea








Magic Mobility +61 3 8791 5600 enquiries@magicmobility.com.au Explore with us magicmobility.com.au











"WE'VE BEEN

CREATING HIGH-QUALITY

WHEELCHAIRS FOR ALMOST 30 YEARS, SO THAT YOU CAN LIVE THE LIFE

YOU CHOOSE."

Ashley Daff R&D Director

We know that small touches make a big difference. Our powerchair range is configurable and modifiable to suit your individual needs.

ROBUST QUALITY

We never compromise on quality. Your powerchair is your mobility, therefore reliability is paramount.



ENGINEERED TO THE LAST BOLT

We listen to and value your feedback. Our team of in-house engineers work side by side with end users and therapists. We're all about innovation; if we think we can make our powerchairs better, we will.



We are proud that our powerchair range is designed and made in our Australian manufacturing facility.