






BEDRAIL MONITOR - FREE STANDING QUICKSTART GUIDE

 English
 Deutsch
 Français
 Español


Distributed by:  **aidacare**
HEALTHCARE EQUIPMENT

INMISA-BEAM[®]
© 2008 Conseng Pty Ltd


- 1** Setting
Paramètre
Setting
Ajuste





- Install **2**
Installer
Anbringen
Instalar




- 3** Connect to power
Relier à la puissance
An die Leitung anschließen
Conectar con la energía




- 4** Monitor **ON** to set beam
Moniteur **DESSUS** pour placer le faisceau
Monitor **AN** zum Einstellen des Lichtstrahls
Monitor **ENCENDIDO** para fijar la viga




- Adjust beam length
Ajuster la longueur du faisceau
Passen Sie die Strahllänge an
Ajustar la longitud del haz



- 6** Monitor **OFF**
Moniteur **AU LOIN**
Monitor **ABSCHALTEN**
Monitor **APAGADO**



- 7** Connect to Nurse Call
Relier à l'appel d'infirmière
An Krankenschwester-Anruf anschließen
Conectar con la llamada de la enfermera



- 8** Monitor **ON** to set alarm
Moniteur **DESSUS** pour placer l'alarme
Monitor **AN** zum Einstellen der Warnung
Monitor **ENCENDIDO** para fijar la alarma

