

Pulpmatic Eco+ Macerator



- The Pulpmatic Eco+ Medical Pulp Macerator features cross-infection control combined with improved environmental performance and ease of use
- Blockage free processing
- Hands-free open, close and start
- Large chamber capacity
- Reduced water consumption
- Automatic sanitation cycle
- Includes programmable audio messaging
- **Water supply:** 15 mm Cold water 0.5 - 10 bar
- **Soil connection:** 50 mm P Trap to 110 mm waste
- **Operation cycle time:** From 85 seconds / 8 ml EcoCleanse
- **Typical capacity:** 3 large or 6 small pulp products



Compliance

ISO 15883/AS 4187:2019

Warranty Information

- 1 Year

Code	Display Name	Total Length	Total Width	Total Height
UPM100130	Panamatic Eco+ Macerator	500 mm	600 mm	1042 mm