

Spreader Bar - Liko - SlingGuard



- Redesigned SlingGuard sling bar, built to enhance patient safety with a patented design unique to the marketplace
- Quick Release Hook model allows easy attachment / detachment to / from the hoist

 Contract approved

Compliance

AS/NZS ISO 10535:2011

Warranty Information

- 1 year



Code	Display Name	Safe Working Load
LIC359320	SlingGuard 450 - Fixed	363kg
LIC359330	SlingGuard 670 - Twin	363kg
LIC359340	SlingGuard 450 - Quick Hook Multi	363kg
LIC359350	SlingGuard 670 - Twin - Quick Hook Multi Spreader Bar	363kg