


Aspire ComfiMotion Bed Collection - Bed Rail



- Improves bed edge awareness and offers multiple grip points for repositioning and bed mobility
- Designed to be used in conjunction with ComfiMotion bed collection - please read variant names for compatibility
- Extendable Version (BEA830205): Ability to extend to prevent falls and contract to allow safe bed egress and ingress (2nd image)

- ✔ Contracted DVA product
- ✔ Contract approved

Compliance

IEC 60601-2-52:2009



| Code | Display Name |
|-----------|--|
| BEA830210 | Aspire ComfiMotion CARE Bed - Bed Rail |
| BEA830212 | Aspire ComfiMotion Activ Care Bed - Bed Rail |
| BEA830205 | Aspire ComfiMotion ACTIV Care Bed Collection - Extendable Bed Rail |