


## Oxford Stature Patient Lifter



- The Stature lift recognises the true needs of the modern care environment, able to cope with even the most demanding and technical of patient handling situations
- Safe working load of 227 kg and great lifting ranges available (min 390 mm, max 1550 mm)
- The Stature can lift patients from the floor as well as accommodate transfers to higher surfaces, while allowing the choice of different sling systems
- Careful attention to product design ensures correct ergonomic usage, protecting both the patient and caregiver from injury
- Manual or powered cradles are available to assist with repositioning of patients
- The powered base allows the caregiver to open and close the legs without bending or twisting
- Ergonomic foot push-pad to assist with initiating movement
- Item Weight: 73.5 kg

### Warranty Information

- 2 Years

 Contract approved

**ndis**

**Item Code:**

05\_123603161\_0105\_1\_2

**Hoist –**

**Composite:**05\_420\_0105\_1\_2

Click [here](#) for more information

**Compliance**

AS/NZS ISO 10535:2011

Code	Display Name	Total Length	Total Width	Total Height	Lift Range	Leg Gait Range	Electric Leg	Safe Working Load	Material
LSS377373	Stature Patient Lifter	1390 mm	1360 mm	1945 - 2130 mm	390 - 1550 mm	640 - 1230 mm	True	227 kg	Aluminium