

## Snowdon Comfort Bed Package



The Snowdon Bedroom selection offers a highly capable assistive bed in a stylish array of interior design options to complement your home. The ultralow minimum height creates a best practice solution for reducing falls risk. This makes the bed particularly helpful for users with Dementia or other cognitive considerations

A range of natural timber palettes through to soft and vibrant upholstered finishes ensures the homely feel and surroundings of the room are maintained while improving safety and independence.

Selecting your perfect bedroom starts with our range of [custom headboards and finishes](#).

Despite the breadth of aesthetic options, each suite naturally addresses the key [clinical considerations](#) important to maintaining user independence and remains compatible with a full range of [accessories](#).

The Aidacare Bedroom Selection is manufactured to the highest standards and the generous warranty terms reflect the market leading [performance and durability](#) of our products

\*Sublime Tweak shown with Shown with laminate side trim

### Warranty Information

- 10 Years Frame
- 3 Years Electrics
- 1 Year Handset & Cables



Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour
BEB048600	Snowdon Comfort Bed Package – Single - Sublime Teak	2130 - 2330 mm	964 mm	100 - 750 mm	275 kg	105 kg	Sublime Teak
BEB050400	Snowdon Comfort Bed Package – Single - Milano Walnut	2130 - 2330 mm	964 mm	100 - 750 mm	275 kg	105 kg	Milano Walnut
BEB048610	Snowdon Comfort Bed Package – King Single - Sublime Tweak	2130 - 2330 mm	1124 mm	100 - 750 mm	275 kg	115 kg	Sublime Teak
BEB050410	Snowdon Comfort Bed Package – King Single - Milano Walnut	2130 - 2330 mm	1124 mm	100 - 750 mm	275 kg	115kg	Milano Walnut