

Phoenix Chair



- The Phoenix is the most supportive and functional of the Seating Matters chairs
- Available in multiple manual & electric combinations
- Designed to give maximum support to the upper body through lateral and head supports, the Phoenix will also adapt to the changing needs of a person being rehabilitated or meet the long term needs of a progressive condition
- The unique adjustability allows the Phoenix to be used with individuals of different dimensions and with varying needs which is ideal for use in acute wards or changing environments
- The adjustable lateral, shoulder and head supports can conform to the person's individual posture and position with fast and easy adjustments
- With tilt in space and adjustable components, the Phoenix will meet a range of individual support needs, sizes, and levels of mobility making it ideal for multiple users
- Used in rehabilitation settings, it can assist in the development and maintenance of independence or for those with progressive conditions; it can provide additional support as necessary
- Effective for adults with the following disabilities and conditions: Cerebral palsy, Spina bifida, Spinal cord injury, Parkinson's disease, Motor neurone disease, Stroke, Multiple sclerosis, Alzheimer's disease, Advanced dementia & Osteoporosis
- Effective for seating clients with the following postural conditions: Pressure ulcers, Reduced upper body tone, Kyphosis, Scoliosis, Muscle contractures & Oedema

Warranty Information

- 7 year frame



Contract approved

ndis

Item Code:

05_180939182_0105_1_2

Click [here](#) for more information

Compliance

ISO13934-1 / ISO13937-1 / BS6807 / BS7175

Accessories

Click on the cart icon to display the range of optional accessories

- 1 year upholstery

Code	Display Name	Seat Depth	Seat Width	Total Height	Backrest Height	Seat To Floor Height	Safe Working Load	Unit Weight
CHP197805	Chair Pressure Relief - Phoenix Chair - Large - Manual - 160kg	390 - 480 mm	400 - 500 mm	1320 mm	850 mm	600 mm	160 kg	67 kg
CHP197810	Chair Pressure Relief - Phoenix Chair - Large - Electric Leg & Tilt - 160kg	390 - 480 mm	400 - 500 mm	1320 mm	850 mm	600 mm	160 kg	67 kg
CHP197815	Chair Pressure Relief - Phoenix Chair - Large - Electric Leg, Recline & Tilt - 160kg	390 - 480 mm	400 - 500 mm	1320 mm	850 mm	600 mm	160 kg	67 kg