

## Aspire Deluxe General Purpose Sling - Mesh



- Australian standards tested and registered with the TGA
- Versatile sling offering body contoured fabric shaping for maximum patient comfort and support
- Quality - manufactured in the UK using European materials, computer guided stitching and CNC fabric cutting
- Durable - loop webbing has a breaking strain in excess of one tonne and is designed so that it will not tangle or “cork-screw” during cleaning
- Serial Numbered - slings are embroidered with a unique serial number to maintain traceability and compliance at all times
- Service & Maintenance - service record labels ensure you can quickly and easily identify service history
- Easy Sling Identification - all slings are colour coded by size and have reinforced edging for longevity
- Easy Clean - machine washable up to 85°C temperature
- For FIM score 1 or 2

- ✓ Contracted DVA product
- ✓ Contract approved

### ndis

#### Item Code:

05\_123603164\_0105\_1\_2

#### Sling – Composite:

05\_421\_0105\_1\_2

Click [here](#) for more information

### Compliance

AS/NZS ISO 10535:2011

## Warranty Information

- 1 Year

| Code      | Display Name                                 | Safe Working Load | Material |
|-----------|----------------------------------------------|-------------------|----------|
| LSS390716 | Deluxe General Purpose Sling - Mesh - Small  | 200 kg            | Mesh     |
| LSS390717 | Deluxe General Purpose Sling - Mesh - Medium | 200 kg            | Mesh     |
| LSS390718 | Deluxe General Purpose Sling - Mesh - Large  | 300 kg            | Mesh     |

| Code      | Display Name                                             | Safe Working Load | Material |
|-----------|----------------------------------------------------------|-------------------|----------|
| LSS390719 | Deluxe General Purpose Sling - Mesh - Extra Large 300 kg |                   | Mesh     |