

## Aspire ComfiMotion Fixed Companion Bed



- Matching single bed without any adjustability designed to sit flush when positioned together with the ComfiMotion CARE Bed and ComfiMotion LUXE Bed
- It facilitates a side-by-side sleeping arrangement at a fixed height of 370 mm, without compromising on the users care and assistance requirements

- ✔ Contracted DVA product
- ✔ Contract approved

### Compliance

IEC 60601-2-52:2009

### Warranty Information

- 1 Year



Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
BEB046800	Aspire ComfiMotion Fixed Companion Bed - Portobello Grey - Long Single (Headboard Sep)	2030 mm	920 mm	370 mm	180 kg