


## DDC Bed Pan Washer Disinfector - Retrofit



- The Data Logger system is suitable for retrofitting to existing assets already installed by Aidacare as well as incorporated into new build sanitisers by the DDC factory
- Chemical Monitoring can be retro fitted
- Please contact us for more information



### Compliance

ISO 15883/AS 4187:2019



Code	Display Name
UPA200300	Retrofit Chemical Monitoring Kit Fits Midi/Optima 2
UPA200310	Retrofit Chemical Monitoring Kit Fits Maxi/Optima 3
UPA200420	Retrofit Kit + Data Logger + Dual Probes Fits Maxi/Optima 3
UPA200430	Retrofit Kit + Data Logger + Dual Probes Fits Midi/Optima 2