

## Oppo 1083 Wrist Wrap with Thumb Loop



- Provides comfortable compressions and support for weak or stressed wrists
- Helps treatment and prevention of injury
- Retains warmth, which may aid and speed up recovery period
- Universal size

✓ Contracted DVA product

✓ Contract approved

### Support @ Home

**Function:** Upper limb orthoses

**Category:** Prescribed

Click [here](#) for more information

| Code      | Display Name  | Wrist Circumference |
|-----------|---|---------------------|
| FBS301368 | Oppo - Wrist Wrap - Thumb Loop - Neoprene - (UNIVERSAL) | 120 - 230 mm        |