

Oxford Advance Folding Patient Lifter



- Mobile patient lifter for transfer to a seated or lying position
- Compact design can be folded easily for storage or transportation
- Epoxy powder coated frame
- Single speed handset control
- Emergency lowering
- Over-sized push handles for excellent manoeuvrability
- Manual foot operated leg gating and tapered leg design for close access
- Suits a large range of patient slings
- Supplied with battery and charger
- Item Weight: 31.7 kg

✓ Contracted DVA product

✓ Contract approved

ndis

Item Code:

05_123603161_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Warranty Information

- 2 Years



Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Total Length	Total Width	Total Height	Lift Range	Leg Gait Range	Safe Working Load	Material
LSS377300	Advance Folding Patient Lifter	1300 mm	550 mm	1360 - 1860 mm	350 - 1860 mm	570 - 1000 mm	155 kg	Aluminium