

DONJOY DP151EB TRIZONE ELBOW SLEEVE



- Three compression zones deliver peak performance.
- Silicone banding injected directly over the sleeve.
- Targets pain from golf and tennis elbow.
- Stabilizes the joint and supports the bicep and triceps muscles.
- All-natural Carbonized Bamboo is thermal regulating and anti-bacterial to help eliminate odours.
- Anti-migration technology helps to keep the elbow sleeve secured to the bicep and eliminates slippage.
- Ergonomic fit.

- ✓ Contracted DVA product
- ✓ Contract approved

| Code | Display Name |
|-----------|---|
| FBS306100 | DONJOY DP151EB TRIZONE ELBOW SLEEVE - SML |
| FBS306090 | DONJOY DP151EB TRIZONE ELBOW SLEEVE - MED |
| FBS306080 | DONJOY DP151EB TRIZONE ELBOW SLEEVE - LGE |
| FBS306110 | DONJOY DP151EB TRIZONE ELBOW SLEEVE - XL |