

Aspire Overbed Table - Thermoform Recessed Top



- Recessed Thermoform top to prevent spills
- Lightweight mobile table positions overbed or chair
- Can be used for eating, drawing or other activities
- Flat top ideal for hospital or home use
- Fully height adjustable
- “C” base shape fitted with four castors enables close positioning to the bedside or chair
- Powder coated mild steel frame construction
- Cream standard laminate top with vinyl wrap around the edge, or “thermoform” cream top which is heat resistant

✔ Contracted DVA product

✔ Contract approved

ndis

Item Code:

05_180315322_0103_1_2

Minor Bed – Composite:

05_428_0103_1_2

Click [here](#) for more information

Warranty Information

- 1 Year



Compliance

ISO134185:2003 / EN ISO
13485:2012

Code	Display Name	Total Length	Total Width	Total Height
BEA028000	Aspire Overbed Table - Thermoform Recessed Top - Cream	430 mm	800 mm	730 - 1095 mm