

## DONJOY AIRCAST AC-05 PNEUMATIC TENNIS ELBOW SUPPORT WITH AIR PAD



- Focused compression for less constriction.
- Breathable material enhances comfort and wearability.
- Simple, effective design makes it easy to apply.
- One size fits most for reduced inventory requirements. Application:
- Tennis elbow
- Golfer's Elbow

✓ Contracted DVA product



Code	Display Name
FBS306070	DONJOY AIRCAST AC-05 PNEUMATIC TENNIS ELBOW SUPPORT WITH AIR PAD - UNIVERSAL