

AGED CARE WAREHOUSE CLEARANCE

ALL PRODUCTS
BETWEEN
50-70% OFF

SALE
HURRY,
WHILST
STOCK
LASTS!

Item Number	Item Description	WAS	NOW	QTY
BEA007489	Alarm - Cura - Soft bed pad - 6mths (760mm x 510mm)	\$ 120.00	\$ 36.00	8
BEA007502	Alarm - Cura - Break Cord USB Mini to RJ10 - for Universal Bed / Chair Pads / Floor Mats	\$ 27.50	\$ 8.25	3
BEA007691	Alarm - Alerta - Standard Cable - Ring/Sleeve	\$ 51.00	\$ 15.30	1
BEA014200	Bed Pole - Grip Handle	\$ 241.00	\$ 72.30	2
BEA014605	Bed Cradle - Aspire - Fixed - Treated Steel ^	\$ 69.00	\$ 20.70	36
BEA024123	Bed End Set - CS7 wide - Amherst 42" w/ Ctrl Panel & T-burg - Solar Oak	\$ 869.00	\$ 260.70	1
BEA024126	Bed Extender - Invacare - 100mm - for Invacare All CS Beds - Adjustable Width	\$ 550.00	\$ 165.00	1
BEA024128	Assist Rail - Invacare - for Invacare CS7 Wide Deck Bed	\$ 402.00	\$ 120.60	1
BEA024129	Assist Bar - Invacare - for Invacare CS7 Bed	\$ 323.00	\$ 96.90	1
BEA024131	Self Help Pole - Invacare - for Invacare All CS Beds	\$ 509.00	\$ 152.70	4
BEA028100	Overbed/Chair Table - Aspire - Wide C Base - Thermoform Recessed Split Tilt Top - Cream	\$ 349.00	\$ 174.50	2
BEA029215	Emergency Battery Back Up - Medicraft	\$ 243.00	\$ 72.90	1
BEA029705	Self Help Pole - Aidacare AC/Ward Beds	\$ 184.00	\$ 92.00	3
BEA029814	Moulded ABS Bed Ends - Aidacare Ward Bed/FL250 ACUTE	\$ 260.00	\$ 78.00	1
BEA030050	Bed Mover/Transporter - Aidacare FL250 Bed (Set of 2))	\$ 809.00	\$ 404.50	3
BEA827090	Bedside Cabinet - Aspire - 1 Drw 1 Dr - Compact - 735 x 530 x 505mm	\$ 975.00	\$ 487.50	1
BEA827120	Bedside Cabinet - Aspire - 1 Drw 1 Dr - Compact - 755 x 500 x 530mm	\$ 550.00	\$ 275.00	1
BEA827150	Bedside Cabinet - Aspire - 3 Drw - Laminate - 615 x 600 x 445mm	\$ 495.00	\$ 247.50	1
BEB045325	Bed - Kerry - Vienna - Single 3 Section w/ 125mm CL Castors	\$ 4,449.00	\$ 1,334.70	5
BEB045535	Bed Invacare 1600 4 section 125mm Castors CL (HB and FB- NOT included)	\$ 2,469.00	\$ 740.70	2
BEM046450	Mattress - Aspire Lifecomfort - Standard - Narrow - 2000x850x130mm - 120kg	\$ 519.00	\$ 259.50	2
BEM046725	Mattress - Aspire Lifecomfort Premier Care - Single - 180kg	\$ 869.00	\$ 434.50	3
BEM046820	Mattress - Lifecomfort with Safety Sides 1940x900x130mm (Cover Only)	\$ 266.00	\$ 133.00	1
BEM047704	Mattress Cover - Protect A Plas - Heavy Duty 1900x910x180mm (7in)	\$ 57.00	\$ 28.50	2
BTC062500	Shower Commode Seat - K Care - H/Made - Closed 550mm	\$ 839.00	\$ 251.70	1
BTC069530	Shower Commode Seat - 4HC - Closed Front w/ Pan Carrier Runners - 550mm Wide	\$ 469.00	\$ 140.70	1
BTC069590	Shower Commode Seat - 4HC Verto TIS - Closed Front - 700mm Wide	\$ 739.00	\$ 221.70	1
BTC071600	Folding Safety Arms - for K Care Shower Commodes	\$ 513.90	\$ 154.17	1
BTC076500	Bowl Carrier - for K Care Junior 400mm Shower Commode - S/Steel	\$ 293.00	\$ 87.90	1
BTC079600	Contoured Headrest - for K Care Shower Commodes	\$ 999.00	\$ 299.70	1
BTC087630	Chest belt	\$ 109.00	\$ 32.70	2
BTC087636	Anti Tippers - Invacare - for Aquatec Ocean & VIP Shower Commodes	\$ 123.00	\$ 36.90	1
BTC092331	Whizard winged Deflector (for Raz AP/SP/CAT/ART/Zum/AP600/SP600/T600)	\$ 124.00	\$ 37.20	1
BTS116100	Shower Chair - Wall Mounted - Fold Away - K Care - AS Large Padded 125kg	\$ 2,358.57	\$ 707.57	1
BTS118750	Mobile Shower Trolley Hose & Fitting	\$ 201.00	\$ 60.30	4
BTS147130	Shower trolley hose kit	\$ 245.00	\$ 73.50	1
BTT153200	Over Toilet Aid - K Care - 445mm - Zinc	\$ 132.00	\$ 39.60	3
BTT156700	Over Toilet Aid - U Frame Adjustable - 510mm (splashguard sep)	\$ 112.90	\$ 33.87	1
BTT157106	Toilet Grab Rail - Throne - Wrap Around - S/Steel	\$ 739.00	\$ 221.70	1
BTT161505	Other Bathroom - Freedom - Slipper Pan	\$ 16.50	\$ 4.95	143
BTT165100	Toilet Surround - Freedom - HBA418HD - 540mm - 120kg	\$ 85.00	\$ 42.50	1
CHP192900	Chair Highback - Comfort Chair Deluxe - 580mm - pc steel	\$ 629.00	\$ 314.50	1
CHP193210	Chair Depth Adj - Todd Chair - High Care	\$ 1,279.00	\$ 639.50	4
CHP197980	Chair Pressure Care - Seating Matters Sorrento - Manual	\$ 6,803.40	\$ 3,401.70	1
CHP198960	Chair Pressure Relief - Aspire COVE Chair - Grey/Andreas Pumice - Medium - 160kg ^	\$ 3,689.00	\$ 1,844.50	1
CHP198965	Chair Pressure Relief - Aspire COVE Chair - Navy/Bronx Steel - Medium - 160kg ^	\$ 3,689.00	\$ 1,844.50	3
CHP198970	Chair Pressure Relief - Aspire COVE Chair - Grey/Andreas Pumice - Large - 160kg ^	\$ 3,689.00	\$ 1,844.50	24
CHP198975	Chair Pressure Relief - Aspire COVE Chair - Navy/Bronx Steel - Large - 160kg	\$ 3,689.00	\$ 1,844.50	2
CHP208400	Chair Footstool - Oval Tube - Maroon Vinyl	\$ 101.00	\$ 30.30	24
CHP208410	Chair Footstool - Oval Tube - Blue Vinyl	\$ 101.00	\$ 30.30	16
CHP208420	Chair Footstool - Oval Tube - Beige Vinyl	\$ 101.00	\$ 30.30	25
CHP213300	Hilite Chair Rear Adj Legs w/ Wheels - Vanilla/Champ	\$ 201.00	\$ 60.30	3
CHP216300	Chair Footrest - Princess Chair/Bed & Deluxe - Air Comfort	\$ 477.00	\$ 143.10	2
CHP216352	Chair Footrest - Air Comfort - Deluxe V2 - Large old style	\$ 331.00	\$ 99.30	2
LAC353230	Baskets to Suit T416N	\$ 25.00	\$ 7.50	4
LAC353810	Trolley - Kerry - Soiled Linen Double w/ Foot Operated Lid	\$ 557.00	\$ 167.10	1
LAC353840	Trolley - Kerry - Soiled Linen Trolley Double w/ Manually Operated Lid	\$ 929.00	\$ 278.70	1
LAC353860	Trolley - Kerry - Soiled Linen Trolley Double S/S w/ Foot Operated Lid	\$ 1,829.00	\$ 548.70	1
LAC354010	Trolley - Kerry - Moist Linen Trolley Small Spring-Loaded Base Only- Galvanized	\$ 719.00	\$ 215.70	3
LIC350197	Ceiling Hoist - Prism Transactive TXE200DPTXY (C450) - Pwr Traverse motor and XY	\$ 4,649.00	\$ 1,394.70	5
LIC353300	Ceiling Hoist - Prism - Track Curved - 15 Degree	\$ 649.00	\$ 194.70	1
LIC353320	Ceiling Hoist - Prism - Track Curved - 45 Degree	\$ 909.00	\$ 272.70	8
LIC353360	Ceiling Hoist - Prism - Bracket End Stop - Long w/Charger	\$ 416.00	\$ 124.80	1
LIC353370	Ceiling Hoist - Prism - Charging Dock for Transactive Xtra	\$ 374.00	\$ 112.20	1
LIC353380	Ceiling Hoist - Prism - Sonell Charger for Transactive Xtra	\$ 535.00	\$ 160.50	1
LIC355901	Brass Air Tube Inserts - Prism For P200	\$ 1.25	\$ 0.38	7
LIC355904	Hub Strap Guide - Prism For P200 - Thick White	\$ 26.50	\$ 7.95	10
LIC355995	CP200 PCB	\$ 1,199.00	\$ 359.70	1

AGED CARE WAREHOUSE CLEARANCE

ALL PRODUCTS
BETWEEN
50-70% OFF

SALE
HURRY,
WHILST
STOCK
LASTS!

Item Number	Item Description	WAS	NOW	QTY
LIC357055	Ceiling Hoist Griffin Track End Stop	\$ 36.10	\$ 10.83	1
LIC357143	C450 Manual Borealis Lifter & Hanger	\$ 4,369.00	\$ 1,310.70	1
LIC357144	C450 Power RTC Borealis Lifter & Hanger	\$ 4,919.00	\$ 1,475.70	1
LIC357193	Ceiling Hoist Safety Reacher Bar 24"	\$ 300.00	\$ 90.00	1
LIC359150	Carriage Hook for Multirall or Likorall	\$ 187.00	\$ 56.10	2
LIC359160	Adapter for Multirall or Likorall	\$ 43.70	\$ 13.11	2
LIC359270	Extension Arm Adapter - Liko Multirall or Likorall - Fits adapter MR/LR	\$ 91.00	\$ 27.30	6
LIC359330	Spreader Bar - Liko SlingGuard 670 twin - Max 363kg (800lbs)	\$ 839.00	\$ 251.70	2
LIC359350	Spreader Bar - Liko SlingGuard 670 twin with Quick-Hook Multi - Max 363kg (800lbs)	\$ 849.00	\$ 254.70	2
LIC359360	Hand control 2 button - Without display - LikoGuard	\$ 709.00	\$ 212.70	1
LIC359370	Hand control 2 button - With display - LikoGuard	\$ 929.00	\$ 278.70	2
LIC359390	Wall Charger - LikoGuard, AU/NZ	\$ 391.00	\$ 117.30	1
LIC359525	Ceiling Hoist - Likorall 250 ES - Natural 250Kg (550 lbs)	\$ 3,859.00	\$ 1,157.70	1
LIC359560	Hand Control - Likorall ES-4MT (for Lifting and Transfer Motor)	\$ 483.00	\$ 144.90	1
LIC360100	Rail Cover, 5 pcs a 2.0m (79 in) - Liko	\$ 89.00	\$ 26.70	1
LIC360900	UltraTwistTM Wide - incl, hand control. Max 500 kg (1100 lbs).	\$ 3,829.00	\$ 1,148.70	1
LIC361000	Curtain Bracket T 300 mm - T-System - Liko	\$ 72.00	\$ 21.60	2
LSS371600	Lifter Sling - Oxford - Quickfit Deluxe Padded Polyester - Medium	\$ 298.00	\$ 149.00	1
LSS371645	Lifter Sling - Oxford Comfort Sling Polyester - Large	\$ 404.00	\$ 202.00	2
LSS371698	Lifter Sling - Oxford - Quickfit Deluxe Poly - Small	\$ 217.00	\$ 108.50	1
LSS371699	Lifter Sling - Oxford - Quickfit Deluxe Poly - Large	\$ 283.00	\$ 141.50	1
LSS374909	Sling - General Purp Head Support - Fab - Small	\$ 327.00	\$ 163.50	1
LSS374925	Sling - General Purpose Paediatric - Medium	\$ 261.00	\$ 130.50	1
LSS375150	Sling - Full Body Mesh Medium	\$ 423.00	\$ 211.50	1
LSS376100	Sling - Supreme Hygiene - Fabric Small	\$ 327.00	\$ 163.50	1
LSS377415	Lifter Yoke to suit Alto Hoist	\$ 205.00	\$ 61.50	1
LSS377769	Lifter Sling - Allegro - Head Support Pivot - Mesh - Large	\$ 444.00	\$ 222.00	1
LSS378800	Lifter Sling - Standing Lifter - Kerry - Large	\$ 297.00	\$ 148.50	1
LSS380300	Lifter Sling - Kerry - GP w/ Neoprene Belt MED (S.M.S.)	\$ 363.00	\$ 181.50	1
LSS380400	Lifter Sling - Kerry - GP w/ Neoprene Belt LGE (S.M.S.)	\$ 363.00	\$ 181.50	1
LSS380700	Lifter Sling - Gen Purp Mesh - Kerry - Large LR406NE3	\$ 297.00	\$ 148.50	1
LSS382400	Lifter Sling - Full Body Polyester w/ Div Leg - Kerry - Medium LR408N2	\$ 376.00	\$ 188.00	1
LSS386202	Lifter - Invacare - 2 Point Spreader Bar - Jasmine	\$ 262.00	\$ 78.60	1
LSS388805	Limb Lifter - Single Patient Use - Wide - 850mm x 250mm - 120kg per Limb (Box of 10)	\$ 247.00	\$ 74.10	1
LSS388810	Pannus Support - SallySling - Single Patient Use - 120kg (Box of 5)	\$ 191.00	\$ 57.30	1
LSS390674	Integrated Weigh Scale - Allegro	\$ 2,679.00	\$ 803.70	2
LSS393410	Wall Mounted Battery Charger, AU/NZ for CBL Battery.	\$ 679.00	\$ 203.70	1
LSS393455	Wheel Set 100/75, (lower clearance) - Fits Viking L & XL	\$ 250.00	\$ 75.00	1
LSS393792	Quick Release Univeral Hook & Q-Link 13 Assembly - Liko	\$ 197.00	\$ 59.10	1
LSS394160	Sling - SafetyVest Mod 93 - Liko - Med - w/4 Loop Clips Fits Sabina II 350 SlingBar	\$ 809.00	\$ 242.70	2
LSS394170	Sling - SafetyVest Mod 93 - Liko - Lge - w/4 Loop Clips Fits Sabina II 350 SlingBar	\$ 849.00	\$ 254.70	2
LSS394180	Sling - ComfortVest Mod 95 - Liko - Sml - w/4 Loop Clips Fits Sabina II Comfort SlingBar.	\$ 779.00	\$ 233.70	1
LSS394210	Sling - ComfortVest Mod 95 - Liko - XL - w/4 Loop Clips Fits Sabina II Comfort SlingBar.	\$ 849.00	\$ 254.70	1
LSS394240	Sling - HygieneVest Highback Mod 55 - Liko - Poly w/ reinforced leg supports - Lge - 200kg	\$ 1,129.00	\$ 338.70	2
LSS394300	Sling (Walking) - MasterVest - Front Close Mod 64 - Liko - Poly - XXL - 200kg	\$ 929.00	\$ 278.70	1
LSS394734	Sling - Liko - ComfortSling Plus 360 - Mesh - XL	\$ 613.00	\$ 306.50	2
LSS394800	Liko Extension Loop, green, 12 cm (5 in)	\$ 17.40	\$ 5.22	1
MWS449610	WChair - Tilt in Space - Aspire Rehab RS - 460mm - 24in Wheels - 180kg - Gunmetal Grey	\$ 2,479.00	\$ 1,239.50	1
PTA493400	HipSaver - SlimFit - High Compliance - Female - XXL	\$ 141.00	\$ 98.70	1
PTA493900	HipSaver Slim Fit Veterans Kit Female Medium	\$ 529.00	\$ 370.30	2
PTA494053	HipSaver - Slim Fit - Pant Only - Female - X Sml	\$ 65.00	\$ 45.50	1
PTA511400	HipSaver NH Veterans Kit Female Medium	\$ 529.00	\$ 370.30	9
PTA511500	HipSaver NH Veterans Kit Female Small	\$ 529.00	\$ 370.30	2
PTA512810	HipSaver - QuickChange - Starter Kit (3-Pants, 1-Pads) - Female - Lge	\$ 269.00	\$ 188.30	2
PTA512949	HipSaver QC Veterans Kit Male Small	\$ 529.00	\$ 370.30	2
PTA513200	HipSaver - SlimFit - Starter Kit (3-Pants, 1-Pads) - Female - Sml	\$ 269.00	\$ 188.30	1
PTA513802	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Female - XXL	\$ 170.00	\$ 119.00	1
PTA513907	HipSaver - Nursing Home - Veterans High Compliance Pack - Male - Sml	\$ 839.00	\$ 587.30	1
PTA521898	Boot - Sheepskin Hi Temp - Short Slipper - C/Toe - X Sml	\$ 215.00	\$ 150.50	2
PTA533319XL	Limb Protection - DermaSaver - Amputee Stump Cover - XL	\$ 98.00	\$ 68.60	1
PTA533327L	Limb Protection - GelBodies - Heel & Ankle Protector - Lge	\$ 118.00	\$ 82.60	1
PTA534102L	Limb/Cast Protection - Limbo - Adult Below Elbow - Waterproof - Lge	\$ 61.00	\$ 42.70	1
PTA534109RL	HipSaver - Track Pant - Regular Length - Lge	\$ 164.00	\$ 114.80	1
PTC545200	Cushion - ROHO Single Compartment HP 13.25x16.75in (7x9cell) Sensor Ready	\$ 1,119.00	\$ 559.50	1
PTC549600	Cushion - Cloud - Otto Bock - 18x16in	\$ 1,279.00	\$ 639.50	2
PTC550600	Cloud Cell Black 100g	\$ 68.00	\$ 34.00	1
PTC550800	Cloud Cell Blue 75g	\$ 68.00	\$ 34.00	2
PTC558600	Cushion - JAY - GS 14x17in w/ Incontinence Cover	\$ 1,119.00	\$ 559.50	4

AGED CARE WAREHOUSE CLEARANCE

ALL PRODUCTS
BETWEEN
50-70% OFF

SALE
HURRY,
WHILST
STOCK
LASTS!

Item Number	Item Description	WAS	NOW	QTY
PTC567900	Cushion - JAY - Care 18x17.5in	\$ 949.00	\$ 474.50	2
PTC568000	Cushion - JAY - Care 16x15.5in	\$ 949.00	\$ 474.50	1
PTC574100	Cushion Cover - JAY - GS 14x17in - Incontinence	\$ 154.00	\$ 46.20	1
PTC576400	Cushion Cover - JAY J2 - 18x20in - Recline Incontinence	\$ 186.00	\$ 55.80	2
PTC576500	Cushion Cover - JAY J2 - 15.5x16in - Recline Air Exch	\$ 126.00	\$ 37.80	2
PTC589500	Cushion - ROHO Quadtro Select LP 20x16.75in (11x9cell)	\$ 1,119.00	\$ 559.50	1
PTC590200	Cushion - ROHO Quadtro Select HP 12x12in (6x6cell)	\$ 1,119.00	\$ 559.50	2
PTC590400	Cushion - ROHO Quadtro Select HP 12x13in (6x7cell)	\$ 1,119.00	\$ 559.50	1
PTC594100	Cushion - ROHO Cont Sel 15x15in (8x8cell)	\$ 1,119.00	\$ 559.50	1
PTC594945	Cushion - ROHO Hybrid Elite 1R 18x20in	\$ 839.00	\$ 419.50	1
PTC594947	Cushion - ROHO Hybrid Elite 1R 22x20in Sensor Ready	\$ 1,239.00	\$ 619.50	1
PTC594950SR	Cushion - ROHO Hybrid Elite 1R 20x18in w/HD cover Sensor Ready	\$ 1,079.00	\$ 539.50	1
PTC595120	Cushion - ROHO Hybrid Elite Cover 16x16in	\$ 197.00	\$ 59.10	4
PTC595125	Cushion - ROHO Hybrid Elite Cover 18x18in	\$ 197.00	\$ 59.10	2
PTC595130	Cushion - ROHO Hybrid Elite Cover 18x20in	\$ 197.00	\$ 59.10	1
PTC595800	Cushion Cover - ROHO Nexus Spirit 16x18in	\$ 179.00	\$ 53.70	2
PTC595900	Cushion Cover - ROHO Nexus Spirit 18x16in	\$ 115.00	\$ 34.50	1
PTC597000	Cushion Cover - ROHO Standard HP 13.25x16.75in	\$ 153.00	\$ 45.90	2
PTC599640	Cushion Cover - Synergy Solution 16 x 18in	\$ 76.00	\$ 38.00	2
PTC608015	Cushion - Promedicare Domino - 10 x 12in w/ 3D cover	\$ 1,619.00	\$ 809.50	1
PTC608040	Cushion - Promedicare Domino - 16 x 16in w/ 3D cover	\$ 1,619.00	\$ 809.50	1
PTC608045	Cushion - Promedicare Domino - 16 x 18in w/ 3D cover	\$ 1,619.00	\$ 809.50	4
PTC608050	Cushion - Promedicare Domino - 18 x 18in w/ 3D cover	\$ 1,619.00	\$ 809.50	1
PTC608056	Cushion - Promedicare Domino - 19 x 18in w/ 3D cover	\$ 1,619.00	\$ 809.50	2
PTC608152	Cushion - Promedicare Tanto - 14 x 15 w/ 3D cover	\$ 1,509.00	\$ 754.50	2
PTC608153	Cushion - Promedicare Tanto - 14 x 16.5in w/ 3D Cover	\$ 1,509.00	\$ 754.50	2
PTC608155	Cushion - Promedicare Tanto - 16 x 16in w/ 3D Cover	\$ 1,509.00	\$ 754.50	2
PTC608160	Cushion - Promedicare Tanto - 16 x 18in w/ 3D Cover	\$ 1,509.00	\$ 754.50	3
PTC608162	Cushion - Promedicare Tanto - 16.5 x 18in w/ 3D Cover	\$ 1,509.00	\$ 754.50	2
PTC608163	Cushion - Promedicare Tanto - 16.5 x 20in w/ 3D Cover	\$ 1,509.00	\$ 754.50	3
PTC608165	Cushion - Promedicare Tanto - 18 x 18in w/ 3D Cover	\$ 1,509.00	\$ 754.50	3
PTC608172	Cushion - Promedicare Tanto - 19in x 18in w/ 3D Cover	\$ 1,509.00	\$ 754.50	2
PTC608280	Cushion Cover - Domino - 3D- 10 x 12in	\$ 294.00	\$ 88.20	2
PTC608282	Cushion Cover - Domino - 3D - 12 x 15in	\$ 294.00	\$ 88.20	2
PTC608289	Cushion Cover - Domino - 3D - 18 x 18in	\$ 294.00	\$ 88.20	2
PTC608291	Cushion Cover - Domino - 3D - 19 x 18in	\$ 294.00	\$ 88.20	1
PTC608323	Cushion Cover - Tanto - 3D - 16 x 18in	\$ 294.00	\$ 88.20	6
PTC608324	Cushion Cover - Tanto - 3D - 16.5 x 20in	\$ 294.00	\$ 88.20	3
PTC609000	Cushion - Vicair - O2 Adjuster 6 Top cover - 18x16inch	\$ 819.00	\$ 409.50	3
PTC622570	Cushion - Pride Stealth Spectrum Cushion - Foam - 18in x 18in	\$ 467.00	\$ 233.50	1
PTM636391	Mattress - Foam - Invacare - Softform Premier - King Single - 1970x1070x150mm	\$ 1,679.00	\$ 839.50	1
PTM636420	Mattress - ROHO - Mattress Overlay Prodigy (3 sections)	\$ 889.00	\$ 444.50	1
PTM636550	Mattress Overlay - Suzric Bubble Mat	\$ 257.00	\$ 128.50	2
PTM636565	Mattress Cover - Concave Positioning Cover - Alt Air - 2000x880x100	\$ 889.00	\$ 266.70	2
PTM636625	AirLo - Thermal Hyperstretch Microfibre Coverlet	\$ 494.00	\$ 148.20	1
PTM638010	Mattress - Aspire HYBRID - Single - 9 Cell - 1980x895x170mm - 180kg (Pump Sep: PTM638030)	\$ 1,723.00	\$ 861.50	2
PTM638020	Mattress - Aspire HYBRID ACUTE - Single - 9 Cell - 1980x895x170mm - 230kg (Pump Sep: PTM638030)	\$ 1,939.00	\$ 969.50	1
PTM639030	Mattress Cover - Aspire Concave Positioning Cover - 2000x900x125mm Active Air 5 Single	\$ 519.00	\$ 155.70	1
PTS640040	Stockings - JOBST for Men - Knee - 15-20mmHg - C/Toe - Sml - Blk	\$ 70.00	\$ 21.00	4
PTS640150	Stockings - JOBST UltraSheer - Knee - 15-20mmHg - C/Toe - Med - Blk 1pair	\$ 70.00	\$ 21.00	2
PTS640315	Stockings - JOBST UltraSheer - Thigh w/ Grip Top - 15-20mmHg - C/Toe - XL - Blk	\$ 97.00	\$ 29.10	1
PTS640640	Stockings - JOBST UltraSheer - Knee - 20-30mmHg - C/Toe - XL - Blk	\$ 103.00	\$ 30.90	1
PTS640761	Stockings - JOBST UltraSheer - Thigh w/ Grip Top - 20-30mmHg - C/Toe - XL - Natural	\$ 144.00	\$ 43.20	2
PTS640910	Stockings - JOBST UltraSheer - Knee - 30-40mmHg - C/Toe - Sml - Natural	\$ 102.00	\$ 30.60	1
PTS641030	Stockings - JOBST UlcerCare - Replacement Liners (pkt 3) - 17mmHg - Med	\$ 52.00	\$ 15.60	1
PTS641359	Stockings - JOBST Relief - Thigh w/ Belt (Chaps) Right - 20-30mmHg - Unisex - O/Toe - Sml	\$ 201.00	\$ 60.30	2
PTS641593	Stockings - TED - Below-Knee - Anti-embolic - Toe Hole - XXXL / Reg - White	\$ 38.30	\$ 11.49	3
PTS642183LM	Compression Sleeve - Venosan 7002 - Sleeve & Hand - Silicone top - 30-40mmHg - Long - Med	\$ 197.00	\$ 59.10	2
PTS642508	Stocking Applicator - Easy Sim-Slide On/Off aid for Open Toe Stockings Small	\$ 60.00	\$ 18.00	1
PTS642700	Socks - Jobst ActiveWear - Knee High - 15-20mmHg - Closed Toe - White - Small	\$ 89.00	\$ 26.70	1
PTS642710	Socks - Jobst for Men Casual - Knee High - 15-20mmHg - Closed Toe - Black - Medium	\$ 89.00	\$ 26.70	2
TSE674400	Slide Sheet - SlippySally - Single Patient Use - 2.0m x 1.45m (Box of 50)	\$ 699.00	\$ 209.70	2
TSE674420	Repositioning Sheet - Sally Sheet - Single Patient Use - Bariatric 1200mm wide - 454kg (Box of 5)	\$ 474.00	\$ 142.20	1
TSE674440	Air Assisted Lateral Transfer Mat Single Patient Use - 860 W x 1980 L (mm) - (Box of 10)	\$ 2,499.00	\$ 749.70	1
TSE674460	Air Assisted Lateral Transfer Mat - 1.7m Replacement Supply Hose	\$ 445.00	\$ 133.50	1
TSE674470	Air Assisted Lateral Transfer Mat - Transport Stand - Single Basket	\$ 1,439.00	\$ 431.70	1
UPA200100	DDC Detergent Basket - Lockable	\$ 403.00	\$ 120.90	3
UPA200210	DDC Detergent Pump - MAXI, Optima 3	\$ 1,239.00	\$ 371.70	9

AGED CARE WAREHOUSE CLEARANCE

ALL PRODUCTS
BETWEEN
50-70% OFF

SALE
HURRY,
WHILST
STOCK
LASTS!

Item Number	Item Description	WAS	NOW	QTY
UPA200240	DDC- Optima 2 / Midi - Utensil tray + sides 2 pc set	\$ 142.00	\$ 42.60	1
UPP100110	Plasticware - Bed Pan Support for (91-01-005))	\$ 57.00	\$ 17.10	10
WAT773600	FlexiCart Options - Folding Shelf with Flat Top	\$ 390.70	\$ 117.21	1
WAT805400	Instrument Trolley - Kerry - 500x500 - 1 Drawer S/S - Epoxy PC Frame	\$ 580.00	\$ 290.00	1
WAT806000	Instrument Trolley - Kerry - 500x500 - 1 Drawer S/S - S/S Frame	\$ 1,469.00	\$ 734.50	1
WAT806800	Dressing Trolley - Kerry - 800x500 - 2 S/S Shelves - Epoxy PC Frame	\$ 629.00	\$ 314.50	1
WAT814500	FlexiCart Options - Basket and mounting bracket	\$ 266.00	\$ 79.80	1
WAT820100	IV Pole PH - Kerry	\$ 155.00	\$ 46.50	1
WAT820120	IV pole - Mobile - 5 star base - Dia 75mm Grey Rubber Castors - Flatpack	\$ 280.00	\$ 84.00	1
WAT900710	Treatment Cart - Aidacare - 650x480x980mm - Blue - 2 Drawers	\$ 1,989.00	\$ 994.50	3
WAT900735	Emergency Cart - Aidacare - 650x480x980mm - Red - 5 Drawers	\$ 3,229.00	\$ 1,614.50	1
WAT905110	Emerg Cart TRIBUTE 90 PKG - 8H DR311 - BRK AWY Lock & RHS Handle - CPR Bd, IV, CO2, Defib, Waste - R	\$ 3,579.00	\$ 1,789.50	1
WAT907190	Cart Accessory - Aidacare/Tribute - Anesthesia Storage Tower Assy (needs 2 x A&E Mount Brackets 1214	\$ 180.00	\$ 54.00	8
WAT907520	Avalo Keyless Pcb	\$ 390.00	\$ 117.00	1