

# EOFY Warehouse Clearance

Code	Item Description	Supplier Code	RRP	NEW Clearance \$	QTY Available / SOH (01/06/2020)
<b>BEDROOM ACCESSORIES</b>					
BEA013800	Bed Pole - Dual Prong for Single Bed	KA524	<del>52.00</del>	\$ 22.50	19
BEA014800	Bed Step - Post Operative - Double	KH32DP	<del>303.00</del>	\$ 136.38	1
BEA016895	Bed IV Pole - Kerry - for Kerry Beds	MC10V	<del>138.00</del>	\$ 20.00	17
<b>BEDSIDE COMMODES</b>					
BEC046020	Bedside Commode - K Care - Deluxe - 500mm 135kg Champagne - Zinc	KA500ZD50V03	<del>923.70</del>	\$ 332.39	4
BEC046060	Bedside Commode w/ Swing Out Bowl - K Care - 440mm Slate - Phosphate	KA500P05	<del>382.90</del>	\$ 172.30	12
BEC046170	Bedside Commode w/ Swing Up Arms - K Care - Champagne - Zinc	KA500ZSV03	<del>518.00</del>	\$ 232.92	3
<b>SHOWER COMMODES</b>					
BTC061600	Shower Commode - K Care - Platform Footrest - 500mm - S/Steel (seat sep)	KA112S50	<del>1,729.50</del>	\$ 750.00	2
BTC061700	Shower Commode - K Care - Platform Footrest - Maxi 550mm 190kg - S/Steel (seat sep)	KA112S55	<del>2,085.10</del>	\$ 750.00	2
BTC061800	Shower Commode - K Care - Platform Footrest - Maxi 600mm 210kg - S/Steel (seat sep)	KA112S60	<del>2,332.60</del>	\$ 750.00	3
BTC062000	Shower Commode - K Care - SA Footrest - 500mm - S/Steel (seat sep)	KA114S50	<del>1,952.30</del>	\$ 750.00	3
BTC062200	Shower Commode - K Care - SA Footrest - Maxi 600mm 210kg - S/Steel (seat sep)	KA114S60	<del>2,491.90</del>	\$ 750.00	2
BTC062300	Shower Commode - K Care - SA Footrest - Maxi 650mm 225kg - S/Steel (seat sep)	KA114S65	<del>2,784.40</del>	\$ 750.00	1
BTC077900	Shower Commode - K Care - Self Prop - SA Footrest - 440mm - S/Steel (seat sep)	KA120SR	<del>1,765.40</del>	\$ 750.00	3
BTC070600	Shower Commode - K Care - Platform Footrest - 440mm - S/Steel (seat sep)	KA112S	<del>1,181.70</del>	\$ 750.00	7
BTC070800	Shower Commode - K Care - Platform Footrest - 440mm - Zinc (seat sep)	KA112Z	<del>890.00</del>	\$ 500.00	2
BTC070900	Shower Commode - K Care - SA Footrest - 440mm - 160kg - S/Steel (seat sep)	KA114S	<del>1,190.50</del>	\$ 500.00	1
BTC071000	Shower Commode - K Care - SA Footrest - 400mm - S/Steel (seat sep)	KA114SJ	<del>1,667.40</del>	\$ 750.00	5
BTC071100	Shower Commode - K Care - Tilt - 440mm - S/Steel (seat sep)	KF11507	<del>4,439.30</del>	\$ 2,500.00	2
BTC071350	Shower Commode - K Care - Tilt - 500mm - S/Steel (seat sep)	KF11512	<del>5,230.60</del>	\$ 2,500.00	1
<b>SHOWER COMMODES ACCESSORIES</b>					
BTC071700	Folding Safety Arms - K Care - for K Care Shower Commodes - Exchange	KO01115	<del>513.90</del>	\$ 125.00	6
BTC072250	Side Panel - K Care - for K Care KA115SRA - with Screws	KA115SRA-SP	<del>234.00</del>	\$ 50.00	26
BTC072270	Shower Commode Arms - Padded Arms (for KA117 Recline) - K Care - Not Fitted	KA117-SPA	<del>47.00</del>	\$ 10.00	1
BTC072500	Shower Commode Seat - K Care - Plastic Seat Closed, Grey	KA11C	<del>67.00</del>	\$ 15.00	9
BTC075300	Backrest - K Care - Curved Padded - STD Width - for K Care Shower Commodes - Fitted	KA11PBF	<del>310.80</del>	\$ 214.44	4
BTC075500	Shower Commode Seat - K Care - Padded Closed - PU Blue 440mm	KA11CPU	<del>138.10</del>	\$ 95.33	4
BTC076600	Bowl Carrier - for K Care Shower Commodes - S/Steel	KA11SBP	<del>133.90</del>	\$ 92.37	14
BTC076950	Footplate - K Care - for K Care Shower Commodes - Grey	KS90026	<del>30.00</del>	\$ 26.15	38
BTC077100	Pan Carrier - for K Care Shower Commodes	KA11SPP	<del>103.90</del>	\$ 71.64	4
BTC077200	Pelvic/Chest Strap - for K Care Shower Commodes - Not Fitted - Requires Rivetting	KA11ST	<del>81.70</del>	\$ 56.37	17
BTC078300	Anti Tip Bars - for K Care Shower Commodes	KA12B	<del>165.60</del>	\$ 114.28	6
BTC079900	Stump Support - Left - for K Care Shower Commodes	KA12LP	<del>468.30</del>	\$ 323.07	2
BTC081604	Commode Seat - K Care - for K Care Bedside Commode - Padded - Slate	50US	<del>66.00</del>	\$ 25.00	10
BTC085000	Commode Seat - 320mm - for K Care Shower Commodes - H/Made Open	KA11HCOP	<del>370.00</del>	\$ 255.33	1
<b>SHOWER CHAIRS &amp; STOOLS</b>					
BTS105600	Shower Chair w/ Arms - K Care - 500mm Seat - Zinc	KA220ZA50	<del>192.00</del>	\$ 132.51	2
BTS106000	Shower Stool w/ Arms & Padded Seat - K Care - 550mm Seat - Zinc	KA222ZA55	<del>778.40</del>	\$ 362.00	8
BTS114700	Shower Chair w/ Arms - K Care - 440mm Seat - S/Steel	KA220SA	<del>137.00</del>	\$ 50.00	2
BTS114800	Shower Chair no Arms - K Care - 440mm Seat - Zinc	KA220Z	<del>126.40</del>	\$ 87.19	15
BTS116600	Shower Stool Padded Seat Replacement - 350mm - K Care	KA22P	<del>31.00</del>	\$ 10.00	5
<b>TOILET AIDS</b>					
BTT153900	Over Toilet Aid Splash Guard - for K Care KA410 - 445mm	KA41SN	<del>15.00</del>	\$ 10.00	58
BTT154700	Over Toilet Aid - K Care - Folding 445mm - Zinc (splashguard sep)	KA430Z	<del>218.70</del>	\$ 99.00	2
<b>DAY CHAIRS</b>					
CHP191400	Chair Lowback - K Care Utility - 550mm Maxi 175kg - Vanilla/Champ	KA586V0355	<del>921.30</del>	\$ 499.00	1
CHP191500	Chair Lowback - K Care Utility - 650mm Maxi 210kg - Vanilla/Champ	KA586V0365	<del>1,168.70</del>	\$ 499.00	1
<b>PRESSURE RELIEF CHAIRS</b>					
CHP198190	Chair Pressure Relief - Air Comfort - Deluxe V2 Bed - Medium	AC59116	<del>3,289.00</del>	\$ 2,500.00	1
CHP198200	Chair Pressure Relief - Air Comfort - Deluxe V2 Bed - Large	AC59111	<del>3,289.00</del>	\$ 2,500.00	12
CHP198390	Chair Pressure Relief - Air Comfort - My Comfort Chair - 180kg	AC59113	<del>1,971.00</del>	\$ 1,500.00	2
CHP201300	Chair Footrest - Padded Large - K Care - Grey/Slate	KA560LG05	<del>146.00</del>	\$ 125.98	1
CHP210300	Chair Leg Rest - Adjustable - Slate Vinyl	KA58LG05	<del>231.30</del>	\$ 159.02	7
CHP210700	Chair Leg Rest - Adjustable - Champagne Vinyl	KA58LV03	<del>218.80</del>	\$ 150.54	1
CHP212800	Hilite Chair Tray Storage Bracket Black	KA58TS	<del>49.70</del>	\$ 5.00	2
CHP216351	Chair Footrest - Air Comfort - Deluxe V2 - Sml/Med/Lge - centre mount	A059055	<del>231.00</del>	\$ 198.83	18
CHP216352	Chair Footrest - Air Comfort - Deluxe V2 - Large old style	A059040	<del>231.00</del>	\$ 198.84	22
CHP216353	Chair Tray - Air Comfort - Deluxe V2 - Large	A059035	<del>223.00</del>	\$ 184.47	15

For all EOFY order inquiries please call 1300 216 898 or email [online@aidacare.com.au](mailto:online@aidacare.com.au)

CHP216400	Chair Tray Small - Princess Chair/Bed - Air Comfort	A059010	<del>149.00</del>	\$ 30.00	15
CHP217100	Chair Footrest - Air Comfort - Tilt Bed	A059008	<del>280.00</del>	\$ 50.00	1
DLG272065	Replacement plastic tray for Kitchen Trolley	KA35TR	<del>34.00</del>	\$ 10.00	3
FPP307200	Tilt Table Acc - Chest Strap	KH33C	<del>463.60</del>	\$ 90.00	1
FPP307400	Tilt Table Acc - Waist Strap	KH33W	<del>463.60</del>	\$ 90.00	2
<b>LAUNDRY &amp; CLEANING</b>					
LAC353430	Cleaners Trolley & Garbag	KH510	<del>1,291.60</del>	\$ 150.00	1
LAC353580	Clean Linen Trolley Front Cover with Zip	KH522FCZ	<del>433.70</del>	\$ 50.00	12
LAC353730	Trolley - Kerry - Soiled Linen Trolley Single Epoxy Coated	KH540	<del>224.50</del>	\$ 100.00	4
LAC353740	Trolley - Kerry - Soiled Linen w/ Foot Operated Lid	KH540L	<del>362.20</del>	\$ 100.00	9
LAC353750	Trolley - Kerry - Soiled Linen & Garbage - Single w/ Foot Operated Lid	KH540LG	<del>451.70</del>	\$ 100.00	11
LAC353770	Trolley - Kerry - Soiled Linen Trolley Single S/S	KH540S	<del>337.90</del>	\$ 100.00	1
LAC353780	Trolley - Kerry - Soiled Linen Trolley Single S/S w/ Foot Operated Lid	KH540SL	<del>769.10</del>	\$ 150.00	6
LAC353810	Trolley - Kerry - Soiled Linen Double w/ Foot Operated Lid	KH541L	<del>557.00</del>	\$ 100.00	15
LAC353820	Trolley - Kerry - Soiled Linen & Garbage - Double w/ Foot Operated Lid	KH541LG	<del>897.80</del>	\$ 100.00	2
LAC353850	Trolley - Kerry - Soiled Linen Trolley Double S/S	KH541S	<del>739.10</del>	\$ 100.00	4
LAC353890	Trolley - Kerry - Soiled Linen Triple w/ Foot Operated Lid	KH542L	<del>1,422.00</del>	\$ 250.00	1
LAC353915	Trolley - Kerry - Soiled Linen Trolley Triple S/S	KH542S	<del>1,175.30</del>	\$ 150.00	5
LAC354060	Trolley - Kerry - Moist Linen Trolley Spring Base (option) - Large	KH548SBC	<del>728.40</del>	\$ 150.00	1
LAC354110	Trolley - Kerry - Soiled Linen Trolley Lid Plastic Green	54LG	<del>40.80</del>	\$ 20.00	56
LAC354130	Trolley - Kerry - Soiled Linen Trolley Lid Plastic Red	54LR	<del>40.80</del>	\$ 20.00	22
LAC354150	Trolley - Kerry - Soiled Linen Trolley Lid Plastic Yellow	54LY	<del>40.80</del>	\$ 20.00	50
LAC354200	Trolley - Kerry - Laundry Trolley Wet/Dry with Spring Base Large	KH552B	<del>2,419.30</del>	\$ 250.00	1
<b>MANUAL HANDLING</b>					
LSS364600	Hoist Legs for Multi Lift KH401, KH404F	404L	<del>261.00</del>	\$ 195.00	1
LSS378300	Lifter - Standing Lifter - Kerry - Quick Stand - 200kg	KH400G	<del>5,800.50</del>	\$ 2,900.00	2
LSS378500	Hoist Leg Strap & Clip On Strap	KH400LS	<del>156.10</del>	\$ 20.00	4
LSS378700	Lifter Sling - Standing Lifter - Kerry - w/ Buckle LR4005FB	KH4005FB	<del>237.00</del>	\$ 150.00	3
LSS381600	Lifter Sling - Gen Purp Polyester & Head Support - Kerry - Medium LR407N2	KH407N2	<del>279.00</del>	\$ 150.00	1
LSS383700	Weigh Scale - Kerry Lifters - SWL 400kg	KH40EWS	<del>2,787.60</del>	\$ 1,500.00	2
WAF704400	Forearm Walker - K Care - Walking Tutor - Child - 600-790mm	KA360C	<del>728.50</del>	\$ 150.00	2
WAF704500	Forearm Walker - K Care - Walking Tutor - Child w/ Glides - 600-790mm	KA360CG	<del>664.40</del>	\$ 99.00	4
WAF708400	Walking Frame Adjustable Folding Swivel Castors (S.M.S.)	KA353SW	<del>196.10</del>	\$ 60.00	1
WAF708900	Forearm Walker - K Care - Walking Tutor - Adult - 850-1350mm	KA360A	<del>626.50</del>	\$ 350.00	11
WAF709150	Forearm Walker - K Care - Walking Tutor - Adult Heavy Duty 175kg - 850-1350mm	KA360AHD	<del>1,017.40</del>	\$ 550.00	2
WAF709200	Forearm Walker - K Care - Walking Tutor - Junior - 670-960mm	KA360J	<del>671.80</del>	\$ 350.00	5
WAF709300	Forearm Walker - K Care - Walking Tutor - Junior w/ Glides - 670-960mm	KA360JG	<del>625.10</del>	\$ 350.00	2
WAS745900	Walking Stick Aluminium Crook Handle Adjustable Silver (S.M.S.)	KA320AI	<del>22.00</del>	\$ 5.00	1
WAS747200	Walking Stick Quad 6"x 8" Base Silver (S.M.S.)	KA332Q68	<del>52.00</del>	\$ 20.00	20
<b>WARD &amp; CLINIC</b>					
WAT773600	FlexiCart Options - Folding Shelf with Flat Top	KH22FSF	<del>390.70</del>	\$ 150.00	4
WAT804300	Instrument Trolley Mayo	KH135	<del>1,717.00</del>	\$ 450.00	1
WAT805300	Instrument Trolley - Kerry - 500x500 - 2 Shelves S/S - Epoxy PC Frame	KH140	<del>596.40</del>	\$ 150.00	2
WAT805400	Instrument Trolley - Kerry - 500x500 - 1 Drawer S/S - Epoxy PC Frame	KH140D	<del>580.00</del>	\$ 450.00	3
WAT807300	Dressing Trolley - Kerry - 800x500 - 1 S/S Drawer - Chrome Frame	KH151D	<del>967.00</del>	\$ 250.00	3
WAT807400	Dressing Trolley - Kerry - 800x500 - 2 S/S Drawers - Chrome Frame	KH151DD	<del>2,002.70</del>	\$ 450.00	1
WAT807800	Dressing Trolley - Kerry - 500x500 - Flat Top - Stainless Steel	KH152	<del>985.60</del>	\$ 150.00	3
WAT809100	Medicine Cart - Kerry - Webster Cart - 3 Drawers (48) - Beige	KH220WB	<del>3,223.80</del>	\$ 800.00	1
WAT814900	FlexiCart Options - Chart Holder - Standard A4	KH22CHL	<del>232.00</del>	\$ 100.00	2
WAT815700	FlexiCart Options - Drive Unit	KH22DU	<del>7,847.70</del>	\$ 1,500.00	1
WAT817950	Flexicart Recessed 30mm Clear Plastic Top	KH22PTC	<del>460.60</del>	\$ 50.00	1
WAT818550	FlexiCart Options - Waste Bin w/mounting bracket	KH22W	<del>265.00</del>	\$ 100.00	3
WAT819100	IV Stand - Mobile - Kerry - Polyurethane Base	KH310	<del>410.70</del>	\$ 100.00	4
WAT819200	IV Stand - Mobile - Kerry - Plastic Base	KH310P	<del>304.60</del>	\$ 100.00	2
WAT822900	Post Op Step Single - Kerry - Grey Large	KH325PLG	<del>102.00</del>	\$ 250.00	1
WAT823900	Food Dispenser Trolley Tarp Cover - Kerry - 16 Tray	KH430TC	<del>973.40</del>	\$ 100.00	4
WAT831200	Paediatric Cot Tray - Kerry	KH942T	<del>485.00</del>	\$ 50.00	1
WAT832200	Neonatal Cot - Kerry - Stainless Steel	KH944S	<del>9,544.60</del>	\$ 1,500.00	1

**For all EOFY order inquiries please call 1300 216 898 or email [online@aidacare.com.au](mailto:online@aidacare.com.au)**