

SIGNATURE DUAL ACTION 2 MOTOR



Back Recline



Leg Recline

Colour range



Mink Warwick
Global Microsuede



Lagoon Warwick
Global Microsuede



Heritage Vinyl



Available in
custom upholstery



- Available in three sizes of single and dual motor models
- Dual action recliner – dual control of leg and recline for infinite positioning
- Depth adjustable seat – a choice of two depth settings in 1 inch (25mm) increments for custom comfort and positioning
- Easy to move – four swivel castors facilitate movement. Provided as standard (seat height increases when fitted)
- Features an all-day battery backup to maximise safety and independence
- Comes with integrated lithium battery backup

Code	Size	Seat Height	Seat Depth	Seat Width	Backrest Height	Overall Height	Overall Width	SWL
Mink Warwick – CHP228200 Lagoon Warwick – CHP228210 Heritage – CHP228220	Small	480mm	460 – 385mm	480mm	580mm	980mm	740mm	
Mink Warwick – CHP228230 Lagoon Warwick – CHP228240 Heritage – CHP228250	Medium	480mm	490 – 515mm	510mm	620mm	1050mm	780mm	155kg
Mink Warwick – CHP228260 Lagoon Warwick – CHP228270 Heritage – CHP228280	Large	500mm	510 – 535mm	540mm	650mm	1110mm	780mm	