

MASS QLD ORDER FORM SUMMARY PAGE

Effective:

ORDERING INFORMATION: Phone 1300 133 120

DEALER INFORMATION

Date	_____	Name	_____
Account Code	_____	Address	_____
Sales Rep	_____	City / State	_____
Mobile	_____	Postcode	_____
		Fax	_____
Signature	_____	Email	_____

CUSTOMER INFORMATION

Customer Name	_____	Contact	_____
MASS Plaque #	_____	Address	_____
MASS PO#	_____	City / State	_____
Signature	_____	Postcode	_____
Rep Name & Ph #	_____	Phone	_____
		Fax	_____
		Email	_____


PRESCRIBER INFORMATION

Prescriber Name	_____	Organisation	_____
Mobile	_____	Address	_____
		City / State	_____
		Postcode	_____
Signature	_____	Phone	_____
		Fax	_____
		Email	_____


PRODUCT INFORMATION

The Aspire EVOKE 2 and EVOKE 2 HD are highly adjustable, lightweight, aluminium wheelchairs designed to enable a personalised fit for every user. A comprehensive range of compatible accessories is available for increased independence and enhanced clinical outcomes.


Max. User Weight	S.W.L 140kg, HD S.W.L - 180kg
Chair Weight	14.8 - 16.2kg
Seat Height - adjustable (floor to top of fabric)	440 - 480mm
Seat Width	350mm 400mm 450mm 500mm 570mm 610mm
Seat Depth - adjustable	400 - 475mm
Backrest Height - adjustable (top of seat to top of backrest handle)	430 - 490mm
Overall Width	560mm 610mm 660mm 710mm 780mm 820mm
Armrest Height (top of seat to top of armrest pad)	205 - 225mm
Approvals	AUST R ARTG 291670 Crash tested to AS/NZS3696.196



SEAT & BACKREST HEIGHT ADJUSTABLE



SEAT DEPTH ADJUSTABLE



BACKREST ANGLE ADJUSTABLE

Effective:

ORDERING INFORMATION: Phone 1300 133 120

ACCOUNT INFORMATION

Invoice to _____	Deliver to _____
Account # _____	Name _____
Address _____	Address _____
_____	_____
Email _____	Email _____
Phone _____	Phone _____
Order Date _____	Marked for _____
PO # _____	Order # _____
Number Required _____	_____

EVOKE STANDARD MASS PACKAGE

Colour: Silver	Armrest: Swing-back & Removable	Seat: Depth & Height Adjustable
Frame: Folding Lightweight Aluminum	Wheels: Quick Release, Solid 24" Tyres, Height & Depth Adjustable	Backrest: Height, Tension, Angle & Depth Adjustable
Legrest: Swing-away & Removable	Castors: 8", Height & Angle Adjustable	Upholstery: Dacron Fire Retardant Fabric
Footrest: Height & Angle Adjustable	Brakes: Easy Access Push-to-Lock Brakes	Cushion: Contoured Overlay is optional, Standard cushion Included

FOR OPTIONAL EXTRAS, PLEASE CHOOSE ONE OF THE PACKAGES BELOW:

CHAIR CUSTOMISATION without cushion - PLEASE CIRCLE BELOW

All customisations will be quoted as labour in addition to the base price, customisations can be carried out by Aidacare or by an authorised Aidacare Network Partner.

Seat Width	350mm	400mm	450mm	500mm	570mm	610mm
-------------------	--------------	--------------	--------------	--------------	--------------	--------------

ACCESSORIES - PLEASE SPECIFY QUANTITY

Code	Elevating Legrests	Price	Qty	Code	Chest Supports	Price	Qty		
WCA848100	Left - 350/400mm	_____	_____	WCA848430	Left	_____	_____		
WCA848110	Right - 350/400mm	_____	_____	WCA848440	Right	_____	_____		
WCA848120	Left - 450/500mm	_____	_____	Code Contoured Padded Overlays					
WCA848130	Right - 450/500mm	_____	_____			Price	Qty		
WCA848140	Left - 570/610mm	_____	_____	WCA848330	350mm	_____	WCA848360	500mm	_____
WCA848150	Right - 570/610mm	_____	_____	WCA848340	400mm	_____	WCA848370	570mm	_____
WCA848160	Left - 350/400mm	_____	_____	WCA848350	450mm	_____	WCA848380	610mm	_____
Code Articulating Elevating Legrests				Code Adjustable Armrests					
WCA848170	Right - 350/400mm	_____	_____	WCA848390	Depth & Height Adjustable - Left	_____	_____		
WCA848180	Left - 450/500mm	_____	_____	WCA848400	Depth & Height Adjustable - Right	_____	_____		
WCA848190	Right - 450/500mm	_____	_____	Code Armrest Pads					
WCA848200	Left - 570/610mm	_____	_____	WCA848410	Left	_____	_____		
WCA848210	Right - 570/610mm	_____	_____	WCA848420	Right	_____	_____		
Code Stump Supports				Code Other					
WCA848220	Left - Silver - 350/400/450mm	_____	_____	WCA848500	Rear Wheels w/Drum Brakes Kit - 24in	_____	_____		
WCA848230	Right - Silver - 350/400/450mm	_____	_____	WCA848270	Anti-Tip Bars	_____	_____		
WCA848240	Left - Silver - 500/570/610mm	_____	_____	WCA848280	IV Pole	_____	_____		
WCA848250	Right - Silver - 500/570/610mm	_____	_____	WCA848300	Safety Lap Belt - Metal Seatbelt Buckle	_____	_____		
WCA848260	Extension Bracket - Silver	_____	_____	WCA848470	Extended Back Cane - Pair - 19/21in	_____	_____		
Code Oxygen Bottle Holder				Code Other					
WCA848290	Oxygen Bottle Holder	_____	_____	WCA848480	Headrest Assembly - Deluxe Style	_____	_____		
Code Brake Extension Levers				Code Other					
WCA848310	Left	_____	_____	WCA848580	Calf Strap	_____	_____		
WCA848320	Right	_____	_____	WCA848490	Camber Kit	_____	_____		

MASS Plaque No.	_____
MASS Funded	_____
Client Contribution	_____
TOTAL PRICE	_____